IMAGO-ARG. IO



WEB INTERFACE FOR MANAGING
MAMMOGRAMS PROCESSED WITH
IMAGO ICE REVEAL

The Product



IMAGO-ARG.io is an intuitive web interface specifically designed for healthcare professionals, enabling the uploading, automatic processing, reporting, and downloading of digital mammograms processed using the groundbreaking Imago ICE Reveal* technology.

*FDA 510K on-going

Data Workflow

- Easy Upload: User-friendly interface allowing secure and rapid image uploads.
- Efficient Management: Clear visualization of processing status.
- Secure Download: Generation of encrypted download links ensuring secure access.

Authentication and Authorization

- Secure multi-factor authentication system powered by Firebase.
- Granular and centralized user administration with role-based permissions.

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IMAGO

Control and Administration

- New user registrations exclusively managed by administrators.
- Mandatory verification and authorization.

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
- Automated email notifications.
- Comprehensive security and traceability of platform access and permissions.

Integration with Google Cloud and Healthcare API

- Robust protection and simplified compliance through integration with Google Cloud.
- Automatic anonymization: Guaranteed privacy and ensured regulatory compliance (GDPR/HIPAA).

Robustness from the First Line of Code


- Automatic Backups: Continuous protection and immediate recovery capabilities.
- Performance and Scalability: Modular, robust, and integrated architecture.
- Continuous Auditing: Constant monitoring and complete traceability.
- End-to-End Encryption: Maximum security during data transfer and storage.

```
dicom_service.py X
dicom-utils > dicom_service.py > ...
1 from flask import Flask, request, send_file, jsonify
2 import pydicom
3 from pydicom.dataset import Dataset
4 from io import BytesIO
5 import json
6 import os
7
8 app = Flask(__name__)
9 app.config['MAX_CONTENT_LENGTH'] = 100 * 1024 * 1024
10
11 def load_custom_tag_config():
12     config_path = os.path.join(os.path.dirname(os.path.dirname(__file__)), 'metadata_files', 'cus
13     try:
14         with open(config_path, 'r') as f:
15             return json.load(f)
16     except Exception as e:
17         print(f"Error loading custom tag config: {e}")
18         return {"customTagConfig": {"tags": []}}
19
20 def apply_custom_tags(ds, doctor_uid):
21     # Usar el tag estándar DICOM para Referring Physician Name
22     ds.add_new(0x00080090, 'PN', doctor_uid)
23     print(f"Added Referring Physician Name tag with value: {doctor_uid}")
24
25     # Verificar que el tag fue agregado
26     if 0x00080090 in ds:
27         print(f"Verified Referring Physician Name tag exists with value: {ds[0x00080090].value}")
28     else:
29         print("Warning: Referring Physician Name tag was not added successfully")
30
31 def extract_metadata(ds):
32     metadata = {}
33     fields = [
34         'SOPInstanceUID', 'SeriesInstanceUID', 'StudyInstanceUID',
35         'PatientID', 'PatientName', 'StudyDate', 'Modality',
36         'ImageLaterality', 'ViewPosition'
37     ]
38     for field in fields:
```

PROBLEMS OUTPUT DEBUG CONSOLE TERMINAL PORTS COMMENTS

```
info: [object Object] {"timestamp": "2025-03-24T19:41:53.936Z"}
info: [object Object] {"timestamp": "2025-03-24T19:41:54.288Z"}
info: [object Object] {"timestamp": "2025-03-24T19:41:54.334Z"}
```

Thanks!